

Evaluating lysozyme deposition on contemporary daily disposable contact lenses in a novel *in vitro* blink model

Vivian WY Chan,¹ Chau-Minh Phan,¹ William Ngo,¹ Lyndon Jones¹
¹Centre for Ocular Research & Education, School of Optometry & Vision Science, University of Waterloo, Canada

Introduction

- Current *in vitro* eye models used in contact lens (CL) deposition studies are still rudimentary.
- Most are performed in a simplistic, static glass vial with the lens submerged in a solution.^{1,2,3}
- Limits the ability to mimic the complexity of the ocular environment.^{1,2,3}
- Our *in vitro* model is designed to physiologically mimic representative tear flow, tear volume, air exposure, and the mechanical rubbing produced during blinking.

Purpose

- To evaluate total lysozyme deposition on daily disposable (DD) CLs using a novel *in vitro* eye-blink model.

Methods

- An *ex vivo* lysozyme deposition study was used as a reference for the blink model⁴ (Table 1).
- Three conventional hydrogel (CH) (etafilcon A, omafilcon A, nelfilcon A) and three silicone hydrogel (SH) (delefilcon A, senofilcon A, somofilcon A) DD CL materials were tested.
- The device blink rate was set to 6 blinks/min with a tear flow rate of 1 μ L/min using an artificial tear solution (ATS) containing lysozyme and other typical tear film components.
- After incubation of 2, 4 or 8 hours the lenses were removed, and the lysozyme extracted using acetonitrile: 0.2% trifluoroacetic acid (n=3).

- A separate experiment was conducted with lenses incubated in a vial containing 480 μ L of ATS on an orbital shaker at 60 rpm for 8 hours (n=3).
- Lysozyme activity was measured using a spectrophotometric assay.




Figure 1 *In vitro* eye blink model used in this study. (A) Polyvinyl alcohol (PVA) eyeball (B) PVA eyelid with 3D-printed support structure (C) 3D-printed lower eyelid (D) tubing connector for fluid input.

Results

- Etafilcon A deposited significantly higher amounts of active lysozyme (402 ± 102 μ g/lens, Figure 2) compared to other lens materials after 8 hours ($p < 0.0001$, Figure 3).
- Etafilcon A had a higher amount of active lysozyme using the blink model (Figure 1) compared to the static vial at 8 hours ($p = 0.0435$).
- Somofilcon A ($p = 0.0076$) and senofilcon A ($p = 0.0019$) had a higher amount of lysozyme activity in the vial compared to the blink model.

Table 1 Active lysozyme deposition comparison between an *ex vivo* study⁴ and the blink model on etafilcon A DD CLs.

Active lysozyme deposition (μ g/lens)	2 hours	4 hours	8 hours
<i>Ex vivo</i> study ⁴	122 ± 1	225 ± 5	342 ± 14
Blink model	137 ± 14	201 ± 155	401 ± 102




Figure 2 Active lysozyme per lens for etafilcon A over 8 hours on the blink model.




Figure 3 Active lysozyme per lens over 8 hours for omafilcon A, nelfilcon A, somofilcon A, senofilcon A, and delefilcon A on the blink model.

Conclusions

- The blink model can be tuned to provide quantitative data that closely mimics *ex vivo* studies and can be used to model deposition of lysozyme on CL materials.

References and Funding

1. Bontempo AR, Rapp J. Lipid Deposits on Hydrophilic and Rigid Gas Permeable Contact Lenses. *CLAO J* 1994;20:242-5.
2. Jones L, Senchyna M, Glasier MA, et al. Lysozyme and Lipid Deposition on Silicone Hydrogel Contact Lens Materials. *Eye Contact Lens* 2003;29:S75-9; discussion S83-4, S192-4.
3. Lorentz H, Heynen M, Trieu D, et al. The Impact of Tear Film Components on *In Vitro* Lipid Uptake. *Optom Vis Sci* 2012;89:856-67.
4. Subbaraman LN, Glasier MA, Varikooty J, et al. Protein Deposition and Clinical Symptoms in Daily Wear of Etafilcon Lenses. *Optom Vis Sci* 2012;89:1450-9.

This study was supported by a Mitacs Accelerate grant and the Canadian Optometric Education Trust Fund (COETF).